

uncapping the syringe FIG1-4 the syringe is placed in the Syringe Guide FIG.1-1 in the direction 6. Surface FIG.1-4a comes into contact with surface FIG.4-10 the syringe FIG.1-4 slides forward until the needle FIG.3-8 punctures the center of the vile head FIG.3 magnified view G-7 and the needle can be seen in the vile. During the puncturing process the Syringe Guide-vile assembly has to be in one of the two conditions as mentioned earlier to prevent the vile from slipping out of the Syringe Guide. After the puncturing process has occurred we have the Syringe Guide, vile and syringe assembly FIG.2. To aspirate the medication from the vile, hold the assembly FIG.2 so that the vile is pointing up. With one hand hold the syringe FIG.1-4 and with the other hand hold the plunger of the syringe (no need to hold the vile) and aspirate. The syringe then can be removed by direction opposite FIG.2-6, capped and ready for use. The vile can also be removed by direction opposite FIG.2-5 and stored for later use.

## CLAIMS

I claim:

1. The Syringe Guide consists of a pipe.
2. Each end of the pipe has a different diameter, large and small.
3. Where the large and small diameters meet in the interior surface of the Syringe Guide there is a wall.
4. The Syringe Guide has a slit running its entire length.
5. The slit makes the Syringe Guide flexible.
6. A medication vile can be placed at one end of the Syringe Guide and a syringe can be placed at the opposite end.
7. The vile and the syringe can dynamically move in the Syringe Guide.

8. The vile and syringe, more specifically move towards each other when in the Syringe Guide.
9. In the puncturing process the vile needs to be supported to prevent its movement.
10. The needle of the syringe punctures the vile directly at the center when both in the Syringe Guide.
11. When aspirating the medication from the vile one does not need to hold the vile.
12. The Syringe Guide makes withdrawing medication from a vile easy.
13. The Syringe Guide makes withdrawing medication from the vile safer.
14. The Syringe Guide makes the area between the needle and the head of the vile more sanitary then the outside environment.
15. The Syringe Guide can be made with different materials.
16. The Syringe Guide can be made in different dimensions.
17. The syringe guide can be made in different constructive ways.